## MULINO A TAGLIENTI SM 50

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## MULINO A TAGLIANTI SM 50 IN CUSTOMER USE

Find out which processes our customers use the SM 50 for.

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The research group uses the SM 50 cutting mill to shred plastic waste, biomass, and electronic scrap for the development of new materials for a sustainable world.

*"We use the SM 50 cutting mill to process plastic waste, biomass, and electronic scrap for the development of new materials that support a more sustainable future. The SM 50 allows us to quickly and easily pre-grind a wide variety of samples, enabling mechanochemical processes that ultimately lead to the creation of new materials in a circular mindset — a true asset to our scientific work."*

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
MULINO A TAGLIENTI SM 50


## FASTER AND MORE FLEXIBLE SAMPLE PREPARATION

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### **Different tasks, different users**


In contract laboratories or academic institutions, many people use the same laboratory grinder for a wide variety of samples. It can happen that feathers are ground today and textiles or completely different sample types tomorrow.

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### **Optimized for different sample tasks in a row**

Healthy nutrition is the basis for a happy life. Our food should combine enjoyment and health. To achieve this, it is important that our food is free from harmful substances and that we have precise knowledge of its composition. The quality of agricultural products is tested in environmental and food laboratories.

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### **Grinding of small sample quantities**

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## ARTICOLI

### CUTTING MILL SM 50

**Cutting Mill SM 50 (please order milling chamber set, rotor, hopper and bottom sieve separately)**

20.720.0001		SM 50	200-230V 50/60Hz, with frequency converter for variable speed
20.720.0003		SM 50	110V 50/60Hz, with frequency converter for variable speed

### MILLING CHAMBER SETS SM 50

**Consisting of milling chamber, double acting cutting bars, door insert and 500 ml glass bottle / receptacle**

22.010.0003	For standard applications: Aluminum milling chamber with stainless steel double acting cutting bars, stainless steel cyclone, stainless steel door insert, and 500 ml glass bottle
22.010.0004	For abrasion and corrosion resistant applications: Stainless steel milling chamber with stainless steel double acting cutting bars, stainless steel cyclone, stainless steel door insert, and 500 ml glass bottle
22.010.0005	For grinding without heavy metal contamination: Aluminum milling chamber with steel 1.1740 double acting cutting bars, stainless steel cyclone, heavy metal free steel door insert, and 500 ml glass bottle
22.010.0006	For grinding without cyclone: Aluminium gravity chamber with stainless steel double acting cutting bars, stainless steel door insert, adapter incl. stainless steel 3 l receptacle

### ROTORS SM 50

22.608.0041	3-disc rotor made of stainless steel with reversable cutting tips made of tungsten carbide
22.608.0042	3-disc rotor made of steel 1.0570 with reversable cutting tips made of tungsten carbide, for grinding without heavy-metal contamination
22.608.0039	Parallel-section rotor made of stainless steel with cutting plates made of stainless steel
22.608.0040	Parallel-section rotor made of steel 1.0570 with cutting plates made of steel 1.1740, for grinding without heavy-metal contamination

### HOPPER SM 50

22.785.0011

EasyInspect hopper made of stainless steel with PE plunger

## BOTTOM SIEVES SM 50

### BOTTOM SIEVES MADE OF STAINLESS STEEL 1.4404

03.647.0574  Fori trapezoidali, 0.25 mm

03.647.0575  Fori trapezoidali, 0.5 mm

03.647.0576  Fori trapezoidali, 0,75 mm

03.647.0577  Fori trapezoidali, 1,00 mm

03.647.0578  Fori quadrati, 2,00 mm

03.647.0579  Fori quadrati, 4,00 mm

03.647.0580  Fori quadrati, 6,00 mm

03.647.0581  Fori quadrati, 8,00 mm

03.647.0582 Fori quadrati, 10,00 mm

22.647.0001 Bottom sieves made of stainless steel complete set of 9 bottom sieves 0.25 - 10 mm

### BOTTOM SIEVES MADE OF STEEL 1.0353 FOR GRINDING WITHOUT HEAVY-METAL CONTAMINATION

03.647.0583 Fori trapezoidali, 0.25 mm

03.647.0584 Fori trapezoidali, 0.5 mm

03.647.0585 Fori trapezoidali, 0,75 mm

03.647.0586 Fori trapezoidali, 1,00 mm

03.647.0587 Fori quadrati, 2,00 mm



03.647.0588 Fori quadrati, 4,00 mm

03.647.0589 Fori quadrati, 6,00 mm

03.647.0590	Fori quadrati, 8,00 mm
03.647.0591	Fori quadrati, 10,00 mm
22.647.0002	Bottom sieves made of steel 1.0353 complete set of 9 bottom sieves 0.25 - 10 mm

## ACCESSORIES SM 50

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22.523.0001		Sample bottle 250 ml (10 pcs.)
22.523.0002		Sample bottle 500 ml (10 pcs.)
02.107.0493		Lid for collecting receptacle cyclone 3l (only for cyclone)
02.011.0031		Collecting receptacle 3 litres, stainless steel (for grinding with or without cyclone)
02.183.0146		Adapter for collecting receptacle 3l and 5l (only for SM 50 gravity outlet)
02.011.0029		Collecting receptacle 5 litres, stainless steel (only for SM 50 gravity outlet)

### MILLING CHAMBERS

22.010.0007	Milling chamber for cyclone made of aluminium with stainless steel double acting cutting bars for standard grinding applications (A)
22.010.0008	Milling chamber for cyclone made of aluminium with steel 1.1740 double acting cutting bars for grinding without heavy metal contamination (A)
22.010.0009	Milling chamber for cyclone made of stainless steel with stainless steel double acting cutting bars for corrosion resistant applications (A)
22.010.0010	Milling chamber for gravity use made of aluminium with stainless steel double acting cutting bars (B)
22.730.0004	Door insert made of heavy metal free steel 1.0570
22.730.0003	Door insert made of stainless steel

### STATIONARY CUTTING BARS FOR MILLING CHAMBERS

22.152.0012	Double acting cutting bars made of stainless steel (3 pcs.) for SM 50
-------------	---

22.152.0013 Double acting cutting bars made of steel 1.1740 (3 pcs.) for grinding without heavy metal contamination

22.608.0043 Adjustment tool for double acting cutting bars

#### ACCESSORIES FOR ROTORS

22.908.0001  Reversible cutting tip incl. screws, tungsten carbide, for 3-disc rotor (10 pcs.)

22.151.0010 Cutting plates made of stainless steel for parallel section rotors (3 pcs.)

22.151.0011 Cutting plates made of steel 1.1740 for parallel section rotors (3 pcs.)

#### CYCLONES

22.935.0040 Cyclone made of stainless steel

22.748.0005 Industrial vacuum cleaner for mills, 230 V, 50/ 60 Hz (other electrical versions and spare filter bags available upon request)

#### ADDITIONAL ITEMS

03.725.0047 PE plunger for EasyInspect hopper

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99.200.0045 IQ/OQ Documentation for SM 50