



HOGE ENERGIE KOGELMOLEN E_{MAX}

the revolution in ultrafine grinding

De E_{max} is een geheel nieuw type kogelmolen voor het malen met hoge energie. De unieke combinatie van hoge wrijvings- en impactkrachten resulteert in extreem fijne deeltjes in de kortst mogelijke tijd.


- | sneller en fijner malen dan alle andere kogelmolens
- | toerental van 2000 min⁻¹ veroorzaakt ultrasnel pulveriseren van het monster
- | water koeling laat continu malen toe, zonder tussentijdse onderbrekingen voor het afkoelen
- | temperatuursbewaking tijdens het malen
- | nauwe korrelgrootte-verdeling dankzij een speciaal maalbekerontwerp, waardoor de menggraad van het monster verbeterd wordt



[Klik om video te bekijken](#)

FASTER - FINER - COOLER - THE MOST POWERFUL BALL MILL

- | Max. speed 2000 rpm
- | Up to 5 mm feed size and 0.08 µm final fineness
- | Two grinding stations for jars of min. 50 ml and max. 125 ml
- | GrindControl to measure temperature and pressure inside the jar.
- | Aeration lids to control the atmosphere inside the jar
- | Temperature monitoring and temperature-controlled grinding, water-cooling of jars
- | Storable SOPs and cycle programs, 4 different jar materials for dry and wet grinding




TOEPASSINGSVOORBEELDEN

legeringen, beenderen, koostofvezel, Katalysatoren, cellulose, cement klinkers, keramiek, chemische producten, kleimineralen, kolen, cokes, beton, vezels, glas, gips, ijzererts, kaolin, kalksteen, metaaloxiden, mineralen, ertsen, papier, pigmenten, plantenmaterialen, polymeren, kwarts, halfedelstenen, slib, slakken, bodems, thee, tabak, afvalmonsters, hout, ...

Bezoek onze toepassingsdatabase om de beste oplossing voor uw monstervoorbereiding of analyse te vinden.

UNIQUE GRINDING JAR GEOMETRY

The High Energy Ball Mill Emax combines high-frequency impact, intensive friction, and controlled circular jar movements to a unique and highly effective size reduction mechanism.



EMAX - FUNCTIE & KENMERKEN
INTUITIVE OPERATION



INSERTING THE GRINDING JAR



CLOSING THE JAR CLAMP




OPERATING THE TOUCHSCREEN

FUNCTIONERINGSPRINCIPE

GRIND SIZES IN THE SUBMICRON RANGE

The High Energy Ball Mill Emax combines high-frequency impact, intensive friction, and controlled circular jar movements to a unique and highly effective size reduction mechanism. The grinding jars have an oval shape and are mounted on two discs respectively which move the jars on a circular course without changing their orientation.


The interplay of jar geometry and movement causes strong friction between the grinding balls, sample material and jar walls as well as a rapid acceleration which lets the balls impact with great force on the sample at the rounded ends of the jars. This significantly improves the mixing of the particles resulting in smaller grind sizes and a narrower particle size distribution than is possible to achieve in ball mills.



[Klik om video te bekijken](#)


VOOR VEILIGE EN EFFICIËNTE MAALPROCESSEN

ACCESSORIES FOR MAXIMUM FLEXIBILITY



GRINDING JARS IN 3 DIFFERENT MATERIALS

Available grinding jar sizes are 50 ml, 80 ml and 125 ml, materials include stainless steel, tungsten carbide and zirconium oxide, ensuring contamination-free sample preparation. Grinding balls are available in sizes from 0.1 mm to 25 mm, depending on the material.



[Klik om video te bekijken](#)

AERATION LID (VIDEO)

RETSCH offers a special aeration lid for the grinding jars designed for applications where a special atmosphere is to be maintained in the ball mill jar.



GRINDCONTROL

The GrindControl measures temperature and pressure inside the jar. The system includes a sensor and transmission unit as well as an analysis software.

TECHNISCHE GEGEVENS

| | |
|---|--|
| Toepassingen | nano vermalings, deeltjesverkleining, homogenisatie, mechanisch legeren, colloïdaal malen, verkleining met hoge energie |
| Toepassingsdomein | biologie, bouwmaterialen, chemie, engineering / electronica, geneesmiddelen / farmaceutica, geologie / metallurgie, glas / keramiek, landbouw, milieu/ recycling |
| Toevoermateriaal | medium-hard, hard, breekbaar, vezelig - droog of nat |
| Maalprincipe | impact, wrijving |
| Korrelgrootte materiaal* | < 5 mm |
| Eindfijnheid* | < 80 nm |
| Batchgrootte / Toevoerhoeveelheid* | max. 2 x 45 ml |
| Snelheid bij 50 Hz (60 Hz) | 300 - 2000 min-1 |
| Koeling | geïntegreerde gecontroleerde water koeling optie: externe koeleenheid |
| Temperatuurscontrole | ja (min en max temperatuur kunnen gedefinieerd worden) |
| Aantal maalstations | 2 |
| Maalbekertypes | met geïntegreerde veiligheidssluiting |
| Materiaal van maalwerktuigen | roestvrij staal, wolframcarbide, zirkoniumoxide |
| Maalbekermaten | 50 ml / 125 ml |
| Instelling maalduur | 00:01:00 to 99:59:59 |
| Interval werking | ja, met optionele omkering van de draairichtingreversal |
| Interval tijd | 00:01:00 to 99:59:59 |
| Pauzetijd | 00:01:00 to 99:59:59 |
| Memoriseerbare standaardprocedures | 10 |
| Interface | USB / LAN (RJ45) |
| Aandrijving | 3-fasen asynchroon motor met frequentie omvormer |
| Aandrijfvermogen | 2600 W |
| Gegevens electriciteit | 200-240 V, 50/60 Hz |
| Aansluiting electriciteit | mono fase |
| Beschermingsklasse | IP 30 |
| Vermogen | ~ 3100W (VA) |
| B x H x D gesloten | 625 x 525 x 645 mm |
| Netto gewicht | ~ 120 kg |
| Standaarden | CE |

*afhankelijk van toegevoerd materiaal en configuratie/instellingen van toestel

www.retsch.nl/emax

BESTELGEGEVENS

HIGH ENERGY BALL MILL EMAX

(please order grinding jars and balls [up to 15 mm] separately)

20.510.0001



Emax, 200–240 V, 50/60 Hz, High energy ball mill with 2 grinding stations

SCREW-LOCK GRINDING JARS EMAX

ROESTVRIJ STAAL

01.462.0305 50 ml

01.462.0313 125 ml



WOLFRAMCARBIDE

01.462.0317 50 ml



ZIRCONIUM OXIDE

01.462.0312 50 ml

01.462.0307 125 ml



AERATION LIDS FOR GRINDING JARS EMAX

incl. o-rings and sintered filter (please order lid and grinding jar separately)

22.107.0638 Aeration lid for grinding jars Emax 50ml

22.107.0640 Aeration lid for grinding jars Emax 125ml

Lid insert for grinding jars Emax

03.474.0258 Aeration lid insert for grinding jars Emax 50ml, rostfreier Stahl

03.107.0570 Aeration lid insert for grinding jars Emax 50ml, Zirkonoxid

03.474.0131 Aeration lid insert for grinding jars Emax 50ml, Wolframcarbide

03.474.0260 Aeration lid insert for grinding jars Emax 125ml, rostfreier Stahl

03.107.0565 Aeration lid insert for grinding jars Emax 125ml, Zirkonoxid

DRUK EN TEMPERATUUR MEETSISTEEM GRINDCONTROL

incl. sensors and transmitter unit, case, opening aid and cleaning accessories for MM 500 control / nano / Emax (please order insert of lid and grinding jar separately)

22.782.0032 GrindControl for MM 500 control/nano/Emax grinding jar 125 ml

GRINDCONTROL LID INSERTS

03.474.0242 GrindControl lid insert for MM 500 control/nano and Emax grinding jar 125 ml, stainless steel

03.474.0245 GrindControl lid insert for MM 500 control/nano and Emax grinding jar 125 ml, zirconium oxide

ACCESSORIES FOR GRINDING JARS EMAX

02.486.0051 Klem om maalbekers te openen

22.186.0007 Sintered filter with O-ring, set of 10 pieces

22.864.0001 Valve set M8x1 for GrindControl and aeration lids



05.114.0057 O-ring for grinding jars 50 ml, 1 piece



05.114.0122 O-ring for grinding jars 125 ml, 1 piece



03.362.0036 Cooling lubricant, 100 ml

99.200.0029 IQ/OQ Documentation for Emax

MAALKOGELS






ROESTVRIJ STAAL

22.455.0010 2 mm Ø, 500 g (ong. 110 ml)



| | | |
|-------------|---|--------------------------------|
| 22.455.0011 |  | 3 mm Ø, 500 g (ong. 120 ml) |
| 22.455.0002 |  | 3 mm Ø, 200 stuks (ong. 6 ml) |
| 22.455.0001 |  | 4 mm Ø, 200 stuks (ong. 14 ml) |
| 22.455.0003 |  | 5 mm Ø, 200 stuks (ong. 25 ml) |
| 05.368.0034 |  | 5 mm Ø |
| 05.368.0035 |  | 7 mm Ø |
| 05.368.0063 |  | 10 mm Ø |
| 05.368.0037 |  | 12 mm Ø |
| 05.368.0109 |  | 15 mm Ø |

WOLFRAMCARBIDE

| | | |
|-------------|---|--------------------------------|
| 22.455.0006 |  | 3 mm Ø, 200 stuks (ong. 6 ml) |
| 22.455.0005 |  | 4 mm Ø, 200 stuks (ong. 14 ml) |
| 22.455.0004 |  | 5 mm Ø, 200 stuks (ong. 25 ml) |
| 05.368.0038 |  | 5 mm Ø |
| 05.368.0039 |  | 7 mm Ø |
| 05.368.0071 |  | 10 mm Ø |

05.368.0041



12 mm Ø

05.368.0110



15 mm Ø

ZIRCONIUM OXIDE

32.368.0005



0.1 mm Ø, 0.5 kg (ong. 135 ml)

32.368.0003



0.5 mm Ø, 0.5 kg (ong. 135 ml)

32.368.0004



1 mm Ø, 0.5 kg (ong. 135 ml)

05.368.0089



2 mm Ø, 0.5 kg (ong. 135 ml)

05.368.0090



3 mm Ø, 0.5 kg (ong. 140 ml)

05.368.0146

7 mm Ø

05.368.0094



10 mm Ø

05.368.0096



12 mm Ø

05.368.0113



15 mm Ø