



ПЛАНЕРАТНА ТОПКОВА МЕЛНИЦА РМ 300

The Planetary Ball Mill PM 300 is a powerful and ergonomic benchtop model with two grinding stations for grinding jar volumes up to 500 ml. This setup allows for processing up to 2 x 220 ml sample material per batch. Thanks to the high maximum speed of 800 rpm, extremely high centrifugal forces result in very high pulverization energy and consequently in short processing times.


The PM 300 is ideally suited for virtually all industries where the quality control process places the highest demands on purity, speed, fineness and reproducibility. Thanks to its enormous energy input of up to 64.4 times the acceleration of gravity, this mill is the perfect choice for tasks in research like mechanochemistry (co-crystal screening, mechanosynthesis, mechanical alloying and mechano-catalysis), or for ultrafine colloidal grinding on a nanometer scale.



[Click to view video](#)

THE POWERFUL, ERGONOMIC PLANETARY BALL MILL

- | Max. speed 800 rpm, large sun wheel
- | Up to 10 mm feed size and 0.1 µm final fineness
- | 2 grinding stations for jars of min. 12 ml and max. 500 ml, jars of 12 – 80 ml can be stacked (two jars each)
- | GrindControl to measure temperature and pressure inside the jar.
- | Aeration lids to control the atmosphere inside the jar
- | Bench top model, touch screen, storable SOPs and cycle programs, 5 different jar materials for dry and wet grinding



FAST & POWERFUL

- | Loss-free size reduction down to the submicron range
- | Wet grinding yields particle sizes in the nanometer range (<100 nm)
- | Variable speed from 50 to 800 rpm, speed ratio 1:-2
- | Milling with up to 64.4 x acceleration of gravity
- | Batch-wise processing with max. 2 x 220 ml sample
- | Possibility to stack jars enables simultaneous processing of 4 samples



SAFE AND EASY HANDLING

- | Easy and safe clamping of the grinding jars thanks to lockable sun wheel
- | The Safety Slider prevents starting the machine without securely clamped jars
- | Perfect stability on the lab bench thanks to FFCS -technology
- | Comfortable parameter setting via Touch display
- | Automatic grinding chamber ventilation with directed air stream for perfect jar cooling
- | Programmable starting time
- | Both grinding jars are freely accessible
- | Ergonomic clamping unit with smooth surfaces



SETTINGS & OPTIONS

- | Dry and wet milling possible
- | Suitable for long-term trials, 99 h max.
- | Programmable breaks e.g. for cooling
- | Automatic calculation of the total process time
- | Direction reversal helps to minimize caking effects
- | Smart service interval notification based on usage



REPRODUCIBILITY

- | Reproducible results due to speed control
- | Indication of load-defined speed, independent of set speed
- | 12 SOPs and 4 cycle programs can be stored
- | Power failure backup ensures storage of remaining processing time



THE BEST ALTERNATIVE TO A RETSCH PLANETARY BALL MILL? A RETSCH MIXER MILL.



Benefit from particularly ergonomic handling while achieving the same finenesses down to the nanometer range.

PLANETARY BALL MILL PM 300

**SAFETY FIRST:
AUTOMATED NOTIFICATION OF CLAMPING FORCE**

Operation of the RETSCH planetary ball mills is particularly safe. They feature a robust Safety Slider which ensures that the mill can only be started after the grinding jar has been securely fixed with a clamping device. The self-acting lock ensures that the jar is seated correctly and securely. This proven solid mechanical system is less failure-prone than electronic solutions - the user has full access to the sample at any time. When the electronic system fails, it is not possible to unlock the jars, for example. A unique safety feature of the PM 300 is an acoustic signal and notification in the display when the clamping unit has been fixed with the required force of 25 Nm. RETSCH offers a convenient clamping aid to facilitate the process. This is particularly helpful when the machine is operated in the upper speed range between 600 and 800 rpm.



PLANETARY BALL MILL PM 300

WET AND NANO-SCALE GRINDING WITH THE PM 300

Wet grinding is used to obtain particle sizes below 5 μm , as small particles tend to get charged on their surfaces and agglomerate, which makes further grinding in dry mode difficult. By adding a liquid or dispersant the particles can be kept separated.

To produce very fine particles of 100 nm or less (nano-scale grinding) by wet grinding, friction rather than impact is required. This is achieved by using a large number of small grinding balls which have a large surface and many friction points. The ideal filling level of the jar should consist of 60 % small grinding balls.

For more details on jar filling, wet grinding and sample recovery watch the video.




[Click to view video](#)

The video shows wet grinding in the Planetary Ball Mill PM 100.

TITANIUM DIOXIDE IN 125 ML JAR

The graphic shows the result of pulverizing titanium dioxide (TiO₂) at 650 rpm and 800 rpm in the PM 300 and the net processing time. With the higher energy input at 800 rpm, the particle size decreases faster. However, the increased warming effects at 800 rpm should also be considered, as they might make longer breaks necessary.




Net processing time of titanium dioxide with 0.1 mm grinding balls in sodium phosphate solution

PLANETARY BALL MILL PM 300

APPLICATIONS IN MECHANOCHEMISTRY

RETSCH Planetary Ball Mills are perfectly suited for processes like mechanical alloying or mechanosynthesis. For most reactions, the 1:2 speed ratio of jar to sun wheel of the models PM 100 and PM 200 is fully adequate, as the ball charge produces enough impact energy. However, greater energy is required for some reactions. Here the PM 400 MA can be used with the increased speed ratio of 1:-2.5 or 1:-3.0.

The PM 300 works with a speed ratio of 1:-2, but in contrast to other models, it reaches up to 64.4 x acceleration of gravity thanks to the maximum speed of 800 rpm and the large sun wheel. Together with the option to use four small, stackable grinding jars sized 12 to 80 ml for small scale operations, or two jars sized up to 500 ml for upscaling purposes, this planetary ball mill is highly suitable for research applications in mechanochemistry.



Achievable acceleration in different planetary ball mills dependent on speed setting

PLANETARY BALL MILL PM 300

EASYFIT GRINDING JARS FOR EXCELLENT RESULTS

The performance and the result of sample preparation are also determined by the choice of the grinding jar and its ball charge. The EasyFit range of jars has been specially designed for extreme working conditions such as long-term trials, even at maximum speed of 800 rpm, wet grinding, high mechanical loads and maximum speeds as well as for mechanical alloying. This line of jars is suitable for all RETSCH planetary ball mills.

The new EasyFit grinding jar series features a structure on the bottom of the 50-500 ml jars called Advanced Anti-Twist (AAT). This ensures that the jars are tightly fixed without the risk of twisting, even at high speed, and that wear and tear is drastically reduced. Secure clamping of the jars is made much easier: to find the correct clamping position, a maximum twist of 60° is required.

The geometry of the EasyFit jars in the 50 ml and 250 ml sizes has been enlarged in diameter and reduced in height compared to the previous "comfort" models. This offers two advantages: better grinding results and interchangeable lids, as there are only three diameter dimensions for the entire grinding jar range.

Diameter categories

- | Diameter 1: 12 ml and 25 ml grinding jars
- | Diameter 2: 50 ml, 80 ml and 125 ml grinding jars
- | Diameter 3: 250 ml and 500 ml grinding jars

- | Available jar sizes: 12 ml / 25 ml / 50 ml / 80 ml / 125 ml / 250 ml / 500 ml
- | Innovative Advanced Anti-Twist (AAT) function ensures secure fit of grinding jars
- | High flexibility thanks to suitability of three lid sizes for all seven jar sizes
- | Pressure-tight and dust-proof O-ring sealing prevents material spillage
- | Jars and balls available in 5 materials: hardened stainless steel, tungsten carbide, agate, sintered aluminium oxide, zirconium oxide
- | Stainless steel protective jacket for agate, sintered aluminium oxide, zirconium oxide and tungsten carbide grinding jars
- | A groove between jar body and lid allows for easy opening of the lid, e. g. with the help of a spatula, if there are underpressure effects inside the jar



ADAPTER FOR SPECIAL APPLICATIONS

With a special adapter, co-crystal screening can be carried out in a planetary ball mill, using disposable vials such as 1.5 ml GC glass vials. The adapter features 24 positions arranged in an outer ring with 16 positions and an inner ring with 8 positions. The outer ring accepts up to 16 vials, allowing for screening up to 64 samples simultaneously when using the Planetary Ball Mill PM 400. The 8 positions of the inner ring are suitable to perform trials with different energy input, e.g. for mechanosynthesis research.



[Click to view video](#)

JARS & LIDS FOR SPECIAL APPLICATIONS

- | For colloidal or wet grinding, the use of a grinding jar with a special closure device is recommended
- | The special closure device is designed for ergonomic handling
- | Aeration lids are designed for working under inert atmosphere, for example if oxygen can influence the grinding process or the mechanosynthesis. The lids allow the introduction of gases like argon or nitrogen into the grinding jar.
- | Optional pressure and temperature measuring system PM GrindControl



GrindControl



Aeration lid



[Click to view video](#)

Video:
Aeration lid

Both the aeration lid and GrindControl can now be equipped with inlays of different materials. Thus, the lid can be used for, e. g. a steel and a zirconium oxide jar by simply exchanging the inlay.

PLANETARY BALL MILL PM 300

RECOMMENDED JAR FILLINGS

To produce optimum grinding results, the jar size should be adapted to the sample amount to be processed. The grinding balls are ideally sized 3 times bigger than the largest sample piece. Following this rule of thumb, the number of grinding balls for each ball size and jar volume is indicated in the table below. To pulverize, for

example, 200 ml of a sample consisting of 7 mm particles, a 500 ml jar and grinding balls sized at least 20 mm or larger are recommended. According to the table, 25 grinding balls are required.

| Grinding jar nominal volume | Sample amount | Max. feed size | Recommended ball charge (pieces) | | | | | |
|-----------------------------|---------------|----------------|----------------------------------|-----------|-----------|---------|---------|---------|
| | | | Ø 5 mm | Ø 7 mm | Ø 10 mm | Ø 15 mm | Ø 20 mm | Ø 30 mm |
| 12 ml | до ≤5 ml | <1 mm | 50 | 15 | 5 | - | - | - |
| 25 ml | до ≤10 ml | <1 mm | 95 – 100 | 25 – 30 | 10 | - | - | - |
| 50 ml | 5 – 20 ml | <3 mm | 200 | 50 – 70 | 20 | 7 | 3 – 4 | - |
| 80 ml | 10 – 35 ml | <4 mm | 250 – 330 | 70 – 120 | 30 – 40 | 12 | 5 | - |
| 125 ml | 15 – 50 ml | <4 mm | 500 | 110 – 180 | 50 – 60 | 18 | 7 | - |
| 250 ml | 25 – 120 ml | <6 mm | 1100 – 1200 | 220 – 350 | 100 – 120 | 35 – 45 | 15 | 5 |
| 500 ml | 75 – 220 ml | <10 mm | 2000 | 440 – 700 | 200 – 230 | 70 | 25 | 8 |

If a pulverization process in a planetary ball mill is successful depends on the machine settings but also on the filling level of the grinding jar. The usable volume of the jars depends on the type of sample material. The number of grinding balls given in the table reflect the minimum amount per jar. A better result is obtained with a larger number of suitable balls, if indicated. In exceptional cases, the number of balls can be reduced by not more than 15%; however, that will result in increased wear of the grinding tools.

PLANETARY BALL MILL PM 300

ТИПИЧНИ МАТЕРИАЛИ ЗА ПРОБА


RETSCH planetary ball mills are perfectly suitable for size reduction of, for example, alloys, bentonite, bones, carbon fibres, catalysts, cellulose, cement clinker, ceramics, charcoal, chemical products, clay minerals, coal, coke, compost, concrete, electronic scrap, fibres, glass, gypsum, hair, hydroxyapatite, iron ore, kaolin, limestone, metal oxides, minerals, ores, paints and lacquers, paper, pigments, plant materials, polymers, quartz, seeds, semi-precious stones, sewage sludge, slag, soils, tissue, tobacco, waste samples, wood, etc.

Very hard, abrasive: Industrial diamonds




11 g sample
50 ml tungsten carbide
grinding jars
4 x 20 mm tungsten
carbide grinding balls
4 min at 400 rpm

Hard-brittle: Sodalite mineral



85 g sample
125 ml zirconium oxide
grinding jars
7 x 20 mm zirconium
oxide grinding balls
12 min at 500 rpm

Nano grinding: Aluminum oxide



100 g sample + 190 ml
sodium-phosphate
solution
500 ml zirconium oxide
grinding jars
1 kg 2 mm zirconium
oxide grinding balls
3:30 min net-grinding
at 650 rpm

*Grinding breaks help to
keep temperature
lower*

Medium-hard, tough: Polyester terephthalate resin




125 g sample
500 ml zirconium oxide
grinding jars
8 x 30 mm zirconium
oxide grinding balls
3 min at 350 rpm

PLANETARY BALL MILL PM 300

ПРИНЦИП НА РАБОТА

Смилащите цилиндри са закрепени ексцентрично на слънчевото колело. Посоката на движение на слънчевото колело е обратна на посоката на движение на смилащите цилиндри и имат предавателно съотношение 1:-2.

Смилащите топки в цилиндърите се движат под въздействието на наслагвани ротационни движения, така наречените Кориолис сили. Взаимодействието на енергиите на смилащите топки и цилиндъра, довеждат до много ефективно намаляване размера на частиците.



[Click to view video](#)

PLANETARY BALL MILL PM 300

TECHNICAL DATA

| | |
|---|---|
| Приложения | pulverizing, mixing, homogenizing, colloidal milling, mechanical alloying, mechanochemistry, co-crystal screening |
| Област на приложение | chemistry, биология, геология / металургия, инженеринг / електроника, медицина / фармация, околната среда / рециклиране, селското стопанство, строителни материали, стъкло / керамика |
| Материали | мек, твърд, крехък, влакнести - сухи или мокри |
| Принцип на смилане | удар, триене |
| Максимален размер на захранващият продукт | < 10 мм |
| Големина след смилане* | < 1 микрона, за колоидално смилане < 0.1 микрона |
| Размер на партидата / количество на пробата* | max. 2 x 220 ml |
| Брой на смилачните станции | 2 |
| Съотношение на скоростта | 1 : -2 |
| Скорост на слънчевото колело | 50 - 800 min ⁻¹ |
| Диаметър на слънчевото колело | 180 мм |
| G-force | 64 g |
| Видове смилачни цилиндри | optional aeration covers, safety closure devices |
| Материал на смилачните части | hardened stainless steel, волфрамов карбид, ахат, спечен алуминиев оксид, циркониев оксид |
| Обем на смилачната камера | 12 ml / 25 ml / 50 ml / 80 ml / 125 ml / 250 ml / 500 ml |
| Stackable grinding jars | 12 ml / 25 ml / 50 ml / 80 ml |
| Adapter for single-use glas vials | 24 x 1.5 ml / 7 x 20 ml |
| Настройване времето за смилане прекъсната работа | digital, 00:00:01 to 99:59:59 |
| Време на прекъсванията | да, с обръщане на посоката |
| Време за пауза | 00:00:01 до 99:59:59 |
| Време за пауза | 00:00:01 до 99:59:59 |
| Запаметяване на настройки | 12 |
| Програми за съхранение на цикъла | 4 |
| Възможност за измерване на вложената енергия | да |
| Интерфейс | USB, RASPI |

| | |
|-------------------------------------|---|
| Задвижване | Трифазен асинхронен мотор с честотен преобразувател |
| Мощност на двигателя | 2,5 kW |
| Данни за електрозахранването | 200-240 V, 50/60 Hz |
| Захранване | еднофазово |
| Защита | IP 20 |
| Енергийна консумация | ~ 3335 VA |
| Размери затворена | 745 x 525 x 580 mm |
| Нетно тегло | ~ 118 kg |
| Стандарти | CE |
| Патенти | да |

в зависимост от захранващият материал и настройките на инструмента

www.retsched.com/bg/pm300

ДАННИ ПОРЪЧКА

PLANETARY BALL MILL PM 300

(please order grinding jars and balls separately)

20.570.0001



PM 300 with 2 grinding
stations, speed ratio 1 :
-2

ACCESSORIES PLANETARY BALL MILLS

| | |
|-------------|--|
| 22.661.0005 | Clamping unit for PM 300 |
| 03.025.0178 | Adapter for stacking grinding jars 50 ml - 80 ml |
| 03.025.0182 | Adapter for the use of grinding jars 12 ml and 25 ml (only for PM 300) |
| 03.486.0062 | Opening aid for clamping unit of planetary ball mills |
| 99.200.0041 | IQ/OQ Documentation for PM 300 |

PRESSURE AND TEMPERATURE MEASURING SYSTEM GRINDCONTROL FOR PLANETARY BALL MILLS

**incl. sensors and transmitter unit, insert of lid, software, case, opening aid and cleaning accessories for PM
(please order grinding jars separately)**

| | |
|-------------|---|
| 22.782.0033 | GrindControl for PM grinding jar EasyFit 50 - 125 ml |
| 22.782.0034 | GrindControl for PM grinding jar EasyFit 250 - 500 ml |

GRINDCONTROL LID INSERTS

| | |
|-------------|---|
| 03.474.0243 | GrindControl lid insert for 50, 80, 125 ml, stainless steel |
| 03.474.0246 | GrindControl lid insert for 50, 80, 125 ml, zirconium oxide |
| 03.474.0244 | GrindControl lid insert for 250 or 500 ml, stainless steel |
| 03.474.0247 | GrindControl lid insert for 250 or 500 ml, zirconium oxide |



ACCESSORIES FOR PM GRINDCONTROL WITH GRINDING JARS EASYFIT

| | | |
|-------------|---|---|
| 05.114.0056 |  | O-ring for 50, 80 or 125 ml |
| 05.114.0054 |  | O-ring for 250 ml - 500 ml grinding jars EasyFit (PM) |
| 03.111.0438 | | Flat gasket for 50 ml, 80 ml or 125 ml |
| 03.111.0439 | | Flat gasket for 250 ml - 500 ml |
| 22.186.0007 | | Sintered filter with O-ring, set of 10 pieces |
| 22.864.0001 |  | Valve set M8x1 for GrindControl and aeration lids |

GRINDING JARS EASYFIT

(grinding jars EasyFit are suitable for all planetary ball mills)

HARDENED STAINLESS STEEL

| | | |
|-------------|---|--------|
| 01.462.0239 |  | 12 ml |
| 01.462.0240 |  | 25 ml |
| 01.462.0516 | | 50 ml |
| 01.462.0517 | | 80 ml |
| 01.462.0518 | | 125 ml |
| 01.462.0519 | | 250 ml |
| 01.462.0520 | | 500 ml |

TUNGSTEN CARBIDE

| | | |
|-------------|--|--------|
| 01.462.0494 | | 50 ml |
| 01.462.0495 | | 80 ml |
| 01.462.0527 | | 125 ml |
| 01.462.0497 | | 250 ml |
| 01.462.0498 | | 500 ml |

AGATE

| | |
|-------------|--------|
| 01.462.0509 | 50 ml |
| 01.462.0511 | 80 ml |
| 01.462.0515 | 125 ml |
| 01.462.0502 | 250 ml |
| 01.462.0506 | 500 ml |



SINTERED ALUMINUM OXIDE

| | |
|-------------|--------|
| 01.462.0507 | 50 ml |
| 01.462.0512 | 125 ml |
| 01.462.0499 | 250 ml |
| 01.462.0503 | 500 ml |

ZIRCONIUM OXIDE

| | |
|-------------|--------|
| 01.462.0508 | 50 ml |
| 01.462.0510 | 80 ml |
| 01.462.0513 | 125 ml |
| 01.462.0500 | 250 ml |
| 01.462.0504 | 500 ml |

ADAPTER FOR GLASS VIALS

| | | |
|-------------|---|--|
| 01.462.0540 |  | Adapter for 24 x 1.5 ml glass vials, stainless, hardened steel |
| 22.749.0009 |  | Glass vial 1.5 ml incl. septum cap, 100 pieces |
| 05.181.0112 | | Replacement pressure spring for adapter for 24 x 1.5 ml glass vials, 1 piece |
| 01.462.0541 | | Adapter for 7 x 20 ml glass vials, stainless, hardened steel |
| 22.749.0010 | | Glass vial 20 ml incl. septum cap, 100 pieces |
| 05.181.0044 | | Replacement pressure spring for adapter for 7 x 20 ml glass vials, 1 piece |

ACCESSORIES FOR GRINDING JARS EASYFIT FOR WET GRINDING,

GRINDING WITH INERT ATMOSPHERE AND MECHANICAL ALLOYING (MA)

AERATION LIDS (INCL. INLAY)

| | |
|-------------|---|
| 22.107.0613 | for grinding jars EasyFit 50 ml - 125 ml, hardened stainless steel |
| 22.107.0616 | for grinding jars EasyFit 50 ml - 125 ml, tungsten carbide |
| 22.107.0617 | for grinding jars EasyFit 50 ml - 125 ml, agate |
| 22.107.0615 | for grinding jars EasyFit 50 ml - 125 ml, zirconium oxide |
| 22.107.0618 | for grinding jars EasyFit 250 ml - 500 ml, hardened stainless steel |
| 22.107.0621 | for grinding jars EasyFit 250 ml - 500 ml, tungsten carbide |
| 22.107.0622 | for grinding jars EasyFit 250 ml - 500 ml, agate |
| 22.107.0620 | for grinding jars EasyFit 250 ml - 500 ml, zirconium oxide |
| 22.107.0619 | for grinding jars EasyFit 250 ml - 500 ml, aluminum oxide |
| 22.864.0001 | Spare valve set for aeration lids M8x1 |



INLAY FOR AERATION LID

| | |
|-------------|---|
| 03.474.0225 | for grinding jars EasyFit 50 ml - 125 ml, hardened stainless steel |
| 03.474.0207 | for grinding jars EasyFit 50 ml - 125 ml, tungsten carbide |
| 03.474.0208 | for grinding jars EasyFit 50 ml - 125 ml, agate |
| 03.474.0206 | for grinding jars EasyFit 50 ml - 125 ml, zirconium oxide |
| 03.474.0226 | for grinding jars EasyFit 250 ml - 500 ml, hardened stainless steel |
| 03.474.0210 | for grinding jars EasyFit 250 ml - 500 ml, tungsten carbide |
| 03.474.0211 | for grinding jars EasyFit 250 ml - 500 ml, agate |
| 03.474.0209 | for grinding jars EasyFit 250 ml - 500 ml, zirconium oxide |
| 03.474.0215 | for grinding jars EasyFit 250 ml - 500 ml, aluminum oxide |

AERATION LIDS FOR GRINDING JARS EASYFIT

INCL. O-RINGS AND SINTERED FILTER (PLEASE ORDER LID INSERT AND GRINDING JAR SEPARATELY)

| | |
|-------------|---|
| 22.107.0636 | Aeration lid for grinding jar EasyFit 50 ml - 125 ml |
| 22.107.0637 | Aeration lid for grinding jar EasyFit 250 ml - 500 ml |

INSERT FOR GRINDING JAR EASYFIT

| | |
|-------------|---|
| 03.474.0261 | Aeration lid insert for grinding jar EasyFit 50, 80 oder 125 ml, stainless steel |
| 03.474.0262 | Aeration lid insert for grinding jar EasyFit 50, 80 oder 125 ml, zirconium oxide |
| 03.474.0263 | Aeration lid insert for grinding jar EasyFit 50, 80 oder 125 ml, tungsten carbide |
| 03.474.0268 | Aeration lid insert for grinding jar EasyFit 50, 80 oder 125 ml, agate |
| 03.474.0264 | Aeration lid insert for grinding jar EasyFit 250 oder 500 ml, stainless steel |
| 03.474.0265 | Aeration lid insert for grinding jar EasyFit 250 oder 500 ml, zirconium oxide |
| 03.474.0266 | Aeration lid insert for grinding jar EasyFit 250 oder 500 ml, tungsten carbide |
| 03.474.0267 | Aeration lid insert for grinding jar EasyFit 250 oder 500 ml, aluminum oxide |
| 03.474.0269 | Aeration lid insert for grinding jar EasyFit 250 oder 500 ml, agate |
| 22.186.0007 | Sintered filter with O-ring, set of 10 pieces |
| 22.864.0001 | Valve set M8x1 for GrindControl and aeration lids |



SAFETY CLOSURE DEVICES

| | |
|-------------|---|
| 22.867.0011 | for grinding jars EasyFit 50 ml - 125 ml |
| 22.867.0012 | for grinding jars EasyFit 250 ml - 500 ml |
| 02.486.0055 | Opening aid for safety closure device |

GASKETS FOR GRINDING JARS EASYFIT

O-RINGS

| | |
|-------------|---|
| 05.114.0086 | O-ring for 12 ml grinding jar EasyFit |
| 05.114.0085 | O-ring for 25 ml grinding jar EasyFit |
| 05.114.0054 | O-ring for 250 ml - 500 ml grinding jars EasyFit |
| 05.114.0056 | O-ring for 50 ml - 125 ml grinding jars EasyFit |
| 05.114.0063 | O-ring for 250 ml - 500 ml grinding jars EasyFit, agate |
| 03.111.0438 | Flat gasket for 50 ml, 80 ml or 125 ml |
| 03.111.0439 | Flat gasket for 250 ml - 500 ml |



GRINDING BALLS

HARDENED STEEL

05.368.0029  5 mm Ø

05.368.0030  7 mm Ø

05.368.0059  10 mm Ø


05.368.0032  12 mm Ø


05.368.0108  15 mm Ø

05.368.0033  20 mm Ø


05.368.0057  30 mm Ø

STAINLESS STEEL

22.455.0010  2 mm Ø, 500 g (approx. 110 ml)

22.455.0011  3 mm Ø, 500 g (approx. 120 ml)

22.455.0002  3 mm Ø, 200 pieces (approx. 6 ml)

22.455.0001  4 mm Ø, 200 pieces (approx. 14 ml)

22.455.0003  5 mm Ø, 200 pieces (approx. 25 ml)

05.368.0034  5 mm Ø

05.368.0035  7 mm Ø


05.368.0063  10 mm Ø

05.368.0037  12 mm Ø

05.368.0109  15 mm Ø


05.368.0062  20 mm Ø

05.368.0105  25 mm Ø

05.368.0061  30 mm Ø

TUNGSTEN CARBIDE

22.455.0006  3 mm Ø, 200 pieces (approx. 6 ml)

22.455.0005  4 mm Ø, 200 pieces (approx. 14 ml)

22.455.0004  5 mm Ø, 200 pieces (approx. 25 ml)

05.368.0038  5 mm Ø

05.368.0039  7 mm Ø

05.368.0071  10 mm Ø

05.368.0041  12 mm Ø

05.368.0110  15 mm Ø

05.368.0070



20 mm Ø

05.368.0069



30 mm Ø

AGATE

05.368.0024



5 mm Ø

05.368.0025



7 mm Ø

05.368.0067



10 mm Ø

05.368.0027



12 mm Ø

05.368.0111



15 mm Ø

05.368.0028



20 mm Ø

05.368.0065



30 mm Ø

SINTERED ALUMINUM OXIDE

05.368.0021



10 mm Ø

05.368.0112



15 mm Ø

05.368.0054



20 mm Ø

05.368.0053




30 mm Ø

05.368.0052




40 mm Ø


ZIRCONIUM OXIDE

32.368.0005  0.1 mm Ø, 0.5 kg (approx. 135 ml)

32.368.0003  0.5 mm Ø, 0.5 kg (approx. 135 ml)

32.368.0004  1 mm Ø, 0.5 kg (approx. 135 ml)


05.368.0089  2 mm Ø, 0.5 kg (approx. 135 ml)


05.368.0090  3 mm Ø, 0.5 kg (approx. 140 ml)

22.455.0007  3 mm Ø, 200 pieces (approx. 6 ml)


22.455.0009  5 mm Ø, 200 pieces (approx. 25 ml)


05.368.0146 7 mm Ø


05.368.0094  10 mm Ø

05.368.0096  12 mm Ø

05.368.0113  15 mm Ø

05.368.0093  20 mm Ø

05.368.0106  25 mm Ø

05.368.0092  30 mm Ø