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The screenshot shows the 'Sieving' configuration screen in the Retsch software. The interface is divided into several sections:

- Navigation:** 'Fast sieving', 'Manual sieving', and 'Methods' tabs are visible at the top.
- Search and Actions:** A search bar and buttons for 'Back', 'Print', 'Refresh', and 'Delete' are located below the tabs.
- Define sample information:** This section includes:
  - Characteristics:** 'Sample material\*' (APAZ), 'Sample preparation\*' (No sample preparation), and 'Density' (100 g/cm³).
  - Source:** 'Container'.
  - Sampling:** 'Sample weight' (3000 g).
- Sample tolerance and check:** A toggle switch for 'Weigh-in tolerance' is set to 'On', with a range from -5.00 % to 5.00 %.
- Off:** 'Weighing assistant' is set to 'Off'.
- Comment:** A text input field for additional notes.

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
### SOFTWARE FEATURES

Graphical User Interface, guided sieving	是
Mode of operation	Fast sieving, manual sieving and methods
Languages	German and English
Parameter units changeable	Size: mm/µm/mesh/inch Mass: g/kg Time: s/min Pressure: mbar/Pa/psi Volume: ml/cl/l/cm <sup>3</sup>
SOPs	Yes, integrated in methods, more than 40
Individual sample labelling	Yes, for each sample, reference to method during data choice and in overview
Parameter setting per test sieve	是
LIMS ready	是
possible export format	.xls, .xlsx, .xps, .csv, .pdf, .png, .jpg, .txt, .rtf, html, .mht
Teamviewer online support	是
Screenshot function	是
Back-up function	是
Automatic storage of results	Yes, in machine, auto report function, PDF document or to printer

[www.retsch.cn/AS200jet-pro-pharma](http://www.retsch.cn/AS200jet-pro-pharma)

## ORDER DATA

### AIR JET SIEVING MACHINES AS 200 JET PRO / JET PHARMA

**Air Jet Sieving Machines for test sieves 203 mm (8") Ø, with automatic vacuum regulation and internal balance, incl. lids for test sieves 1" and 2" height and rubber mallet - please order test sieve and MANDATORY vacuum cleaner (Item No. 22.748.0020 or Item No. 22.748.0021) separately!**

30.034.0001		AS 200 100–240 V, 50/60 Hz jet pro	incl. test report acc. to EN 10204 2.2 evaluation software installed (no EasySieve required)
30.034.0002		AS 200 100–240 V, 50/60 Hz jet pharma (no EasySieve required)	incl. test report acc. to EN 10204 2.2 evaluation software installed incl. AuditTrail & User Management


### VACCUUM CLEANER AS 200 JET PRO / JET PHARMA

22.748.0020	Industrial vacuum cleaner ATTIX 33 230V 50/60Hz (other electrical versions available) MANDATORY vacuum cleaner for 30.034.0001 & 30.034.0002 NOT SUITABLE for 32.935.0008
32.524.0017	HEPA filter for ATTIX 33
32.524.0018	Cover for HEPA filter ATTIX 33
32.524.0019	Replacement filter bags for ATTIX 33, 5 pcs.

### CYCLONE AS 200 JET PRO / JET PHARMA

22.935.0042	High-performance cyclone including holder, collecting receptacle (250 ml and 500 ml sample bottle) and 1.5 m FDA certified hose with clamps (industrial vacuum cleaner 22.748.0020 or 22.748.0021 recommended)
05.707.0208	FDA certified hose, 1.5 m
02.183.0151	Connector for cyclone

### ACCESSORIES AS 200 JET / JET PRO / JET PHARMA

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22.523.0001

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22.523.0002

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## TEST SIEVES WITH ELECTROFORMED SHEET

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## ADDITIONAL ITEMS AIR JET SIEVING MACHINES

99.200.0047

IQ/OQ Documentation for AS 200 jet pro

99.200.0046

IQ/OQ Documentation for AS 200 jet pharma

## RECALIBRATION AIR JET SIEVING MACHINES

**The Retsch Air Jet Sieving Machines can be calibrated at Retsch/Germany or at selected subsidiaries/agents**

SER-97.001.0021

Calibration of sieving machines AS 200 jet pro & AS 200 jet pharma  
at Retsch/Germany, incl. test report acc. to EN 10204 2.2